

Guided Surgery

Drill Usage Chart

Description	Part Number	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20
2.0mm Pilot Drill, 17mm Stop	122-017																				
2.0mm Pilot Drill, 21mm Stop	122-021																				
2.0mm Pilot Drill, 24mm Stop	122-024																				
2.0mm Pilot Drill, 28mm Stop	122-028																				
2.5mm Width Increasing Drill, 17mm Stop	CGS-2517																				
2.5mm Width Increasing Drill, 21mm Stop	CGS-2521																				
2.5mm Width Increasing Drill, 24mm Stop	CGS-2524																				
2.5mm Width Increasing Drill, 28mm Stop	CGS-2528																				
2.8mm Width Increasing Drill, 17mm Stop	CGS-2817																				
2.8mm Width Increasing Drill, 21mm Stop	CGS-2821																				
2.8mm Width Increasing Drill, 24mm Stop	CGS-2824																				
2.8mm Width Increasing Drill, 28mm Stop	CGS-2828																				
3.2mm Width Increasing Drill, 17mm Stop	CGS-3217																				
3.2mm Width Increasing Drill, 21mm Stop	CGS-3221																				
3.2mm Width Increasing Drill, 24mm Stop	CGS-3224																				
3.2mm Width Increasing Drill, 28mm Stop	CGS-3228																				
3.7mm Width Increasing Drill, 17mm Stop	CGS-3717																				
3.7mm Width Increasing Drill, 21mm Stop	CGS-3721																				
3.7mm Width Increasing Drill, 24mm Stop	CGS-3724																				
3.7mm Width Increasing Drill, 28mm Stop	CGS-3728																				
4.1mm Width Increasing Drill, 17mm Stop	CGS-4117																				
4.1mm Width Increasing Drill, 21mm Stop	CGS-4121																				
4.1mm Width Increasing Drill, 24mm Stop	CGS-4124																				
4.1mm Width Increasing Drill, 28mm Stop	CGS-4128																				
4.7mm Width Increasing Drill, 17mm Stop	CGS-4717																				
4.7mm Width Increasing Drill, 21mm Stop	CGS-4721																				
4.7mm Width Increasing Drill, 24mm Stop	CGS-4724																				
4.7mm Width Increasing Drill, 28mm Stop	CGS-4728																				
5.4mm Width Increasing Drill, 17mm Stop	CGS-5417																				
5.4mm Width Increasing Drill, 21mm Stop	CGS-5421																				
5.4mm Width Increasing Drill, 24mm Stop	CGS-5424																				
5.4mm Width Increasing Drill, 28mm Stop	CGS-5428																				

Important Note:
All drills should be replaced after they have been used 12 to 20 times, depending upon bone density.

Shop on-line at www.biohorizons.com